



Calibration of Liquid Admixture Report

Date 2-3-26

Customer Name: Daytona Redi mix

Location: Shelby twp plant 3

State: MI

Yard/Plant ID: Plant 3
Vince hagen

DOT#

Product Name	Pulse Rate oz	Control Panel Reading oz.	Calibrated Container Reading oz.	Variance (+/-)%	Pass/Fail
<u>SA-56</u>	<u>15</u>	<u>101</u>	<u>100</u>		<u>P</u>
<u>400mc</u>	<u>1</u>	<u>104</u>	<u>106</u>		<u>P</u>
<u>Dynamon SX</u>	<u>1</u>	<u>99</u>	<u>102</u>		<u>P</u>
<u>Poly R</u>	<u>1</u>	<u>104</u>	<u>104</u>		<u>P</u>
<u>Easy-75</u>	<u>1</u>	<u>100</u>	<u>100</u>		<u>P</u>
<u>V3K</u>	<u>1</u>	<u>101</u>	<u>98</u>		<u>P</u>
<u>Superset +</u>	<u>4.0</u>	<u>328</u>	<u>328</u>		<u>P</u>
<u>Chloride</u>	<u>80</u>	<u>320</u>	<u>320</u>		<u>P</u>

All calibrated container readings are accurate to within 3% of the total amount required.

All air entrainment calibrated container readings are accurate to within 1% of total amount required.

Jordan Rhodes
Jordan Rhodes

DOT or Company Representative Signature: [Signature]